



NRSN2 Antibody

Product Code	CSB-PA016094GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZP1
Immunogen	Human NRSN2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	neurensin 2;NRSN2;C20orf98;FLJ23329;dJ1103G7.6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NRSN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.