



NUDT19 Antibody

Product Code	CSB-PA016168GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8MXV4
Immunogen	Human NUDT19
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	nudix (nucleoside diphosphate linked moiety X)-type motif 19;NUDT19;RP2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NUDT19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.