



NUDT6 Antibody

Product Code	CSB-PA016176GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53370
Immunogen	Human NUDT6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nudix (nucleoside diphosphate linked moiety X)-type motif 6;NUDT6;ASF2GF2;FGF-2;FGF-AS;FGF2AS;bFGF;gfg;gfg-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NUDT6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.