



NUFIP1 Antibody

Product Code	CSB-PA016182GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHK0
Immunogen	Human NUFIP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	nuclear fragile X mental retardation protein interacting protein 1;NUFIP1;NUFIP;bA540M5.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NUFIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.