



NUMB Antibody

Product Code	CSB-PA016186GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49757
Immunogen	Human NUMB
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	numb homolog (Drosophila);NUMB;S171 ;
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NUMB
Target Details	<p>This protein plays a role in the determination of cell fates during development. The encoded protein, whose degradation is induced in a proteasome-dependent manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Four transcript variants encoding different isoforms have been found for this gene.</p>