



OGDH Antibody

Product Code	CSB-PA016307GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02218
Immunogen	Human OGDH
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	oxoglutarate (alpha-ketoglutarate) dehydrogenase (lipoamide); OGDH; AKGDH; E1k; OGDC ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OGDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.