



OXCT2 Antibody

Product Code	CSB-PA017307GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYC2
Immunogen	Human OXCT2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	3-oxoacid CoA transferase 2;OXCT2;FKSG25;FLJ00030;SCOT-T ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OXCT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.