



PAK1IP1 Antibody

Product Code	CSB-PA017405GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWT1
Immunogen	Human PAK1IP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	PAK1 interacting protein 1;PAK1IP1;FLJ20624;MAK11;PIP1;RP11-421M1.5;WDR84;bA421M1.5;hPIP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PAK1IP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.