



# PARD3 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA017444GA01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q8TEW0  |
| <b>Immunogen</b>           | Human PARD3   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human,Mouse,Rat   |
| <b>Tested Applications</b> | ELISA,WB,IHC  |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.  |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Alias</b>               | par-3 partitioning defective 3 homolog (C. elegans);PARD3;ASIP;Baz;Bazooka;FLJ21015;PAR3;PAR3alpha;PARD3A;SE2-5L16;SE2-5LT1;SE2-5T2 ; |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | PARD3   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |