



PARP11 Antibody

Product Code	CSB-PA017459GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR21
Immunogen	Human PARP11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	poly (ADP-ribose) polymerase family, member 11;PARP11;C12orf6;DKFZp779H0122 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PARP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.