



# PARP9 Antibody

<b>Product Code</b>	CSB-PA017472GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8IXQ6
<b>Immunogen</b>	Human PARP9
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	poly (ADP-ribose) polymerase family, member 9; PARP9; BAL; BAL1; DKFZp666B0810; DKFZp686M15238; FLJ26637; FLJ35310; FLJ41418; FLJ43593; MGC:7868 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PARP9
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.