



PCSK4 Antibody

Product Code	CSB-PA017643GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UW60
Immunogen	Human PCSK4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	proprotein convertase subtilisin/kexin type 4;PCSK4;DKFZp434B217;MGC34749;PC4;SPC5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PCSK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.