



PDS5B Antibody

Product Code	CSB-PA017743GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTI5
Immunogen	Human PDS5B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	PDS5, regulator of cohesion maintenance, homolog B (S. cerevisiae);PDS5B;APRIN;AS3;CG008;FLJ23236;KIAA0979;RP1-267P19.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDS5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.