



PEBP1 Antibody

Product Code	CSB-PA017766GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30086
Immunogen	Human PEBP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	phosphatidylethanolamine binding protein 1;PEBP1;HCNP;PBP;PEBP;RKIP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PEBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.