



ECI2 Antibody

Product Code	CSB-PA017768GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75521
Immunogen	Human PECl
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	peroxisomal D3,D2-enoyl-CoA isomerase;PECl;ACBD2;DRS1;HCA88;KIAA0536;dJ1013A10.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ECI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.