



PFN4 Antibody

Product Code	CSB-PA017827GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NHR9
Immunogen	Human PFN4
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	129aa,14kd;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PFN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.