



PGF Antibody

Product Code	CSB-PA017854GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49763
Immunogen	Human PGF
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	placental growth factor;PGF;D12S1900;PGFL;PLGF;PIGF-2;SHGC-10760 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PGF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.