



PGRMC1 Antibody

Product Code	CSB-PA017876GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00264
Immunogen	Human PGRMC1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	progesterone receptor membrane component 1;PGRMC1;HPR6.6;MPR ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PGRMC1
Target Details	Progesterone binding protein is a putative steroid membrane receptor. The protein is expressed predominantly in the liver and kidney.