



PHF11 Antibody

Product Code	CSB-PA017900GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIL8
Immunogen	Human PHF11
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	PHD finger protein 11;PHF11;APY;BCAP;IGEL;IGER;IGHER;NYREN-34;NYREN34 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHF11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.