



PHIP Antibody

Product Code	CSB-PA017921GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWQ0
Immunogen	Human PHIP
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pleckstrin homology domain interacting protein; PHIP; FLJ20705; FLJ45918; MGC90216; WDR11; ndrp ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.