



PHPT1 Antibody

Product Code	CSB-PA017942GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX4
Immunogen	Human PHPT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	phosphohistidine phosphatase 1;PHPT1;CGI-202;DKFZp564M173;HSPC141;PHP14;bA216L13.10 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHPT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.