



PI4KA Antibody

Product Code	CSB-PA017955GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42356
Immunogen	Human PI4KA
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	phosphatidylinositol 4-kinase, catalytic, alpha;PI4KA;FLJ16556;PI4K-ALPHA;PIK4CA;pi4K230 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PI4KA
Target Details	<p>This gene encodes a phosphatidylinositol (PI) 4-kinase which catalyzes the first committed step in the biosynthesis of phosphatidylinositol 4,5-bisphosphate. The mammalian PI 4-kinases have been classified into two types, II and III, based on their molecular mass, and modulation by detergent and adenosine. This protein is a type III enzyme that is not inhibited by adenosine. Two transcript variants encoding different isoforms have been described for this gene.</p>