



PIK3CB Antibody

Product Code	CSB-PA017999GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42338
Immunogen	Human PIK3CB
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phosphoinositide-3-kinase, catalytic, beta polypeptide;PIK3CB;DKFZp779K1237;MGC133043;PI3K;PI3KCB;PI3Kbeta;PIK3C1;p110-BETA ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3CB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.