



PIK3R3 Antibody

Product Code	CSB-PA018006GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92569
Immunogen	Human PIK3R3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phosphoinositide-3-kinase, regulatory subunit 3 (gamma);PIK3R3;DKFZp686P05226;FLJ41892;p55;p55-GAMMA ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3R3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.