



PIKFYVE Antibody

Product Code	CSB-PA018029GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2I7
Immunogen	Human PIP5K3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	phosphatidylinositol-3-phosphate/phosphatidylinositol 5-kinase, type III;PIP5K3;CFD;FAB1;KIAA0981;MGC40423;PIKFYVE;PIP5K ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIKFYVE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.