



PKN2 Antibody

Product Code	CSB-PA018075GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16513
Immunogen	Human PKN2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	protein kinase N2;PKN2;MGC150606;MGC71074;PAK2;PRK2;PRKCL2;PRO2042;Pak-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PKN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.