



PLAGL2 Antibody

Product Code	CSB-PA018119GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPG8
Immunogen	Human PLAGL2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	pleiomorphic adenoma gene-like 2;PLAGL2;FLJ23283 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLAGL2
Target Details	Pleiomorphic adenoma gene-like 2 is a zinc-finger protein that recognizes DNA and/or RNA.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.