



PLEKHG5 Antibody

Product Code	CSB-PA018174GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94827
Immunogen	Human PLEKHG5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pleckstrin homology domain containing, family G (with RhoGef domain) member 5;PLEKHG5;DSMA4;GEF720;KIAA0720;RP4-650H14.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLEKHG5
Target Details	This gene encodes a protein that activates the nuclear factor kappa B (NFKB1) signaling pathway. Mutations in this gene have been found in a family with distal spinal muscular atrophy.