



# PLEKHG5 Antibody

<b>Product Code</b>	CSB-PA018174GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O94827
<b>Immunogen</b>	Human PLEKHG5
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	pleckstrin homology domain containing, family G (with RhoGef domain) member 5;PLEKHG5;DSMA4;GEF720;KIAA0720;RP4-650H14.3 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PLEKHG5
<b>Target Details</b>	This gene encodes a protein that activates the nuclear factor kappa B (NFkB1) signaling pathway. Mutations in this gene have been found in a family with distal spinal muscular atrophy.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.