



PM20D2 Antibody

Product Code	CSB-PA018228GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYS1
Immunogen	Human PM20D2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	peptidase M20 domain containing 2;PM20D2;ACY1L2;bA63L7.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PM20D2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.