



PNCK Antibody

Product Code	CSB-PA018259GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P2M8
Immunogen	Human PNCK
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	pregnancy up-regulated non-ubiquitously expressed CaM kinase;PNCK;BSTK3;CaMK1b;FLJ50403;FLJ50549;FLJ56451;FLJ59811;MGC 45419 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PNCK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.