



PNRC1 Antibody

Product Code	CSB-PA018289GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12796
Immunogen	Human PNRC1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	proline-rich nuclear receptor coactivator 1;PNRC1;B4-2;PNAS-145;PROL2;PRR2;RP11-63L7.5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PNRC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.