



POLR1E Antibody

Product Code	CSB-PA018326GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZS1
Immunogen	Human POLR1E
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	polymerase (RNA) I polypeptide E, 53kDa;POLR1E;FLJ13390;FLJ13970;PAF53;PRAF1;RP11-405L18.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLR1E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.