



# POLR3A Antibody

<b>Product Code</b>	CSB-PA018343GA01HU
<b>Abbreviation</b>	POLR3A
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O14802
<b>Immunogen</b>	Human POLR3A
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	polymerase (RNA) III (DNA directed) polypeptide A, 155kDa;POLR3A;RPC1;RPC155;hRPC155 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	POLR3A