



PPBP Antibody

Product Code	CSB-PA018429GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02775
Immunogen	Human PPBP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	pro-platelet basic protein (chemokine (C-X-C motif) ligand 7);PPBP;B-TG1;Beta-TG;CTAP-III;CTAP3;CTAPIII;CXCL7;LA-PF4;LDGF;MDGF;NAP-2;PBP;SCYB7;TC1;TC2;TGB;TGB1;THBGB;THBGB1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PPBP
Target Details	This protein is a platelet-derived growth factor that belongs to the CXC chemokine family. This growth factor is a potent chemoattractant and activator of neutrophils. It has been shown to stimulate various cellular processes including DNA synthesis, mitosis, glycolysis, intracellular cAMP accumulation, prostaglandin E2 secretion, and sythesis of hyaluronic acid and sulfated glycosaminoglycan. It also stimulates the formation and secretion of plasminogen activator by synovial cells.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.