



PPIL6 Antibody

Product Code	CSB-PA018487GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXY8
Immunogen	Human PPIL6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	peptidylprolyl isomerase (cyclophilin)-like 6;PPIL6;MGC27054;MGC27056;MGC41939;bA425D10.6;dJ919F19.1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPIL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.