



PPP3CA Antibody

Product Code	CSB-PA018572GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08209
Immunogen	Human PPP3CA
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	protein phosphatase 3 (formerly 2B), catalytic subunit, alpha isoform;PPP3CA;CALN;CALNA;CALNA1;CCN1;CNA1;PPP2B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP3CA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.