



# Recombinant *Pseudomonas putida* Phosphatidylserine decarboxylase proenzyme (psd)

<b>Product Code</b>	CSB-EP018033PZN-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A5W9U4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Pseudomonas putida</i> (strain ATCC 700007 / DSM 6899 / BCRC 17059 / F1)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSRLFIIISQ YLLPHHLLSR LAGCIAECRV RWFKNAFTAW FAKRYQVNMS EALVEDLSAY EHFNAFFTRA LKPGARPLDE TPGAILCPAD GAVSQLGPIE HGRIFQAKGH GFSAQELLGG DPAMAAPFMG GEFATIYLSP KDYHRVHMPL AGTLREMVYV PGRLFVSNQT TAENVPELFA RNERVVCLFD TERGPMVAVL VGAMIVASIE TVWAGLVTPP KRELKTFRYD EASRTPIHLE KGAE LGRFKL G
<b>Source</b>	E.coli
<b>Target Names</b>	psd
<b>Protein Names</b>	Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
<b>Expression Region</b>	1-251
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.