



PRPF19 Antibody

Product Code	CSB-PA018761GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMS4
Immunogen	Human PRPF19
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	PRP19/PSO4 pre-mRNA processing factor 19 homolog (S. cerevisiae);PRPF19;NMP200;PRP19;PSO4;SNEV;UBOX4;hPSO4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRPF19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.