



PRUNE2 Antibody

Product Code	CSB-PA018832GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUY3
Immunogen	Human PRUNE2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	prune homolog 2 (Drosophila);PRUNE2;A214N16.3;BMCC1;BNIPXL;C9orf65;DKFZp762K117;FLJ50060;FLJ54876;FLJ59118;KIAA0367;RP11-58J3.2;bA214N16.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRUNE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.