



PSG11 Antibody

Product Code	CSB-PA018851GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ72
Immunogen	Human PSG11
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pregnancy specific beta-1-glycoprotein 11;PSG11;MGC22484;PSG13;PSG14 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PSG11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.