



PSMD13 Antibody

Product Code	CSB-PA018903GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNM6
Immunogen	Human PSMD13
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	proteasome (prosome, macropain) 26S subunit, non-ATPase, 13; PSMD13; HSPC027; Rpn9; S11; p40.5 ;
Immunogen Species	Homo sapiens (Human)
Target Names	PSMD13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.