



PTBP2 Antibody

Product Code	CSB-PA018950GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKA9
Immunogen	Human PTBP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	polypyrimidine tract binding protein 2;PTBP2;FLJ34897;PTB;PTBLP;brPTB;nPTB;nPTB5;nPTB6;nPTB7;nPTB8 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTBP2
Target Details	This protein binds to the intronic cluster of RNA regulatory elements, downstream control sequence (DCS). It is implicated in controlling the assembly of other splicing-regulatory proteins. This protein is very similar to the polypyrimidine tract binding protein but it is expressed primarily in the brain.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.