



PTPMT1 Antibody

Product Code	CSB-PA019023GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUK0
Immunogen	Human PTPMT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	protein tyrosine phosphatase, mitochondrial 1;PTPMT1;DUSP23;FLJ46081;MOSP;PLIP;PNAS-129 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTPMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.