



PXN Antibody

Product Code	CSB-PA019112GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49023
Immunogen	Human PXN
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Protein A Purified
Isotype	IgG
Alias	paxillin;PXN;FLJ16691 ;Paxillin
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PXN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.