



RAB11FIP5 Antibody

Product Code	CSB-PA019159GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXF6
Immunogen	Human RAB11FIP5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RAB11 family interacting protein 5 (class I);RAB11FIP5;DKFZp434H018;GAF1;KIAA0857;RIP11;pp75 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB11FIP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.