



RAB33A Antibody

Product Code	CSB-PA019185GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14088
Immunogen	Human RAB33A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.-20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RAB33A, member RAS oncogene family;RAB33A;MGC1488;RabS10 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB33A
Target Details	This protein belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport.