



RAB3GAP1 Antibody

Product Code	CSB-PA019198GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15042
Immunogen	Human RAB3GAP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RAB3 GTPase activating protein subunit 1 (catalytic);RAB3GAP1;DKFZp434A012;KIAA0066;P130;RAB3GAP;RAB3GAP1 30;WARBM1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB3GAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.