



RAB3IP Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA019201GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96QF0 |
| Immunogen | Human RAB3IP |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | ELISA, WB, IHC |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | RAB3A interacting protein (rabin3); RAB3IP; FLJ14660; FLJ22548; MGC71495; RABIN3 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RAB3IP |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |