



RAB4A Antibody

Product Code	CSB-PA019211GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20338
Immunogen	Human RAB4A
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:2000, IF:1:50-1:500
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	24kd;
Immunogen Species	Homo sapiens (Human)
Target Names	RAB4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.